

# Earn a **\$500** credit towards the purchase of a Focus™ intraoral X-ray when you purchase a Snapshot or Digora® Optime digital intraoral system

Offer valid for purchases made between January 1, 2010 and March 31, 2010  
Offer can not be combined with other promotional offers

## Snapshot - Intraoral Sensor

The Snapshot sensor's innovative, compact electronics combined with CMOS sensor technology make using the Snapshot quick, easy, and flexible.

- ✓ Simple USB connection makes the Snapshot plug and play.
- ✓ Take almost every intraoral image with just one exposure setting. Snapshot's wide dynamic range provides superior images of different parts of the dentition without altering the settings on your X-ray device
- ✓ The new, sophisticated Snapshot sensor design produces up to 20 lp/mm true resolution (26,3 lp/mm theoretical). You get razor-sharp images, with the low noise typical of CMOS technology.

# OR

## DIGORA® Optime

Continuing the SOREDEX heritage of easy-to-use imaging equipment, the DIGORA® Optime takes usability to the next level. The system has self-guiding tutorials with clear visual aids. Even the most novice users will be able to use the system with ease.

- ✓ Touchless operation - The system automatically detects your approaching hand and instantly prepares itself for plate scanning
- ✓ Lightning fast image readout - A high quality size 2 image is read in 5 seconds!
- ✓ Automatic plate detection and erasing - The system automatically detects the size of the imaging plate – no cumbersome cassettes, plate guides or drums are needed. Imaging plates are automatically erased without user interaction and are immediately ready for reuse.
- ✓ DIGORA® Optime - Just like film



For more information contact  
Instrumentarium Dental or Soredex at  
1-800-558-6120

